

FOR CRITICAL PROCESS CLEANING

CANI

COMMITMENT TO CONSTANT AND NEVER ENDING IMPROVEMENT

Complete III One-Step Acid Cleaner And Descaler

For Use On:

Stainless steel
Carbon & urine scales
Oils & proteins
Inorganic salts & scales
Fecal matter
Animal cages
Pharmaceutical/biotech
process cleaning
& derouging

DESCRIPTION

Complete III is a phosphoric acid cleaner/descaler formulated for the research, pharmaceutical, and biotech industries. It effectively removes a wide range of process and research residues, such as inorganic salts and scales, particulate carbon and urine scales, fecal matter, oil and proteins. This product is an effective, ONE-STEP cleaner and descaler in animal cage wash applications. This product is very effective in pharmaceutical/biotech/cosmetic process cleaning and derouging.

Low foaming **COMPLETE III** will not cavitate recirculation pumps. **COMPLETE III** is extremely free rinsing to achieve a clean, residue free surface and has excellent hard water tolerance to eliminate mineral deposits. The surfactant system in **COMPLETE III** remains effective across a wide temperature and pH range, even during pH neutralization prior to discharge, thus reducing redeposition of process residues.

DIRECTIONS

Typical operating concentrations of **COMPLETE III** run from 0.2% to 5.0% (v/v) and depend on temperatures, soil loads, and cleaning methods used. **COMPLETE III** is designed for use in immersion, pressure spray and recirculating spraywash applications and is ideal for animal cage and glassware washing machines and recirculating Clean-In-Place (CIP) systems. Higher concentrations are typically used for periodic derouging and descaling. Contact CANI, Inc. for technical services or product recommendations.



PRODUCT SPECIFICATIONS

Form: Liquid
Odor: characteristic
pH1% solution: 1.8-2.2
Specific Gravity: 1.20-1.24
Solubility: complete
Foam: low

CANI, INC.

1758 Allentown Road, Lansdale, PA 19446

Phone: (610) 222-4500

Fax: (610) 222-4203

CRITICAL PROCESS CLEANING

