

FOR CRITICAL PROCESS CLEANING

CANI

COMMITMENT TO CONSTANT AND NEVER ENDING IMPROVEMENT

Millennium Citrus LF Low Foaming Citrus Cleaner/Degreaser

For Use On:

- Cosmetic Devices
- Pharmaceutical Devices
- Medical Devices
- Food & Beverage processing area

DESCRIPTION

A highly concentrated liquid, low foam citrus cleaner & degreaser designed to cut through tough soil conditions, and leave surfaces sparkling clean with a fresh pleasant citrus odor. Can be used as detergent additive in immersion, pressure washing and circulation (C.I.P.) cleaning for the cosmetic, pharmaceutical, medical device and food and beverage industries.

DIRECTIONS

Low Foam Citrus Cleaner & Degreaser is designed for use in immersion, pressure washing and circulation (C.I.P.) cleaning applications. Typical operating concentrations of **Low Foam Citrus Cleaner & Degreaser** run from 0.2% to 5.0% (v/v), with compatible alkaline cleaners and depend on temperatures, soil loads, and cleaning methods used. Proper care and protective equipment are required for manual cleaning see the MSDS.



PRODUCT SPECIFICATIONS

| | | | |
|----------------------------|-------------|---------------------------------|----------------------|
| Active Ingredients: | 13.0 ± 0.5% | Solubility: | Complete with water |
| Color: | Amber | Emulsification: | Excellent |
| Odor: | Citrus | Detergency: | Excellent |
| pH: | 12.5 - 13.0 | Rinsability: | Complete with water |
| Viscosity: | Water thin | Flash Point: | None (waterbased) |
| Foam: | Low | Biodegradability: | Complete |
| Abrasives: | None | Storage/Stability: | 2 years |
| Phosphate Content: | None | Weight Per Gallon: | 8.75 lbs. per gallon |
| | | Freeze / Thaw Stability: | Keep from freezing |

CANI, INC.

1758 Allentown Road, Lansdale, PA 19446

Phone: (610) 222-4500

Fax: (610) 222-4203

CRITICAL PROCESS CLEANING

